






# 21

Sterilization  
Sterilisation  
Esterilización  
Stérilisation  
Sterilizzazione

Type	Packing		Requirements according to:					
	Description	Symbolic Delivery and Packaging System	DIN EN 285	DIN 58946-1				
			Gravity Procedure	Pre-Vacuum Procedure	Fractional Vacuum Procedure	Fractional Flow Procedure	Steam Injection Procedure	
	Textiles in sterilization container with filter inside lid and bottom.		+	+	+	+	+	+
			+	0	0	+	+	+
	Textiles in sterilization container with filter inside lid, and not perforated bottom.		+	-	-	+	+	+
			+	-	-	+	+	0
	Textiles in sterilization container with filter inside bottom and not perforated lid.		+	0	+	+	+	+
			+	-	0	+	+	+
	Instruments in wire baskets inside sterilisation container with filter inside lid and bottom.		+	+2)	+	+	+	+
	Instruments in wire baskets in sterilization container with filter inside lid and not perforated bottom.		+	-	+2)	+2)	+2)	+3)
	Instruments in wire baskets inside sterilization container with filter inside bottom.		+	0	+	+	+	+

Type	Packing		Requirements according to:					
	Description	Symbolic Delivery and Packaging System	DIN EN 285	DIN 58946-1				
			Gravity Procedure	Pre-Vacuum Procedure	Fractional Vacuum Procedure	Fractional Flow Procedure	Steam Injection Procedure	
	Textiles in sterilization container with valve inside lid and bottom.		<b>+</b>	<b>-</b>	<b>0</b>	<b>+</b>	<b>+</b>	<b>0</b>
			<b>+</b>	<b>-</b>	<b>-</b>	<b>+</b>	<b>+</b>	<b>0</b>
	Textiles in sterilization container with valve inside lid and no valve in the bottom.		<b>+</b>	<b>-</b>	<b>-</b>	<b>+</b>	<b>+</b>	<b>-</b>
			<b>+</b>	<b>-</b>	<b>-</b>	<b>+</b>	<b>+</b>	<b>-</b>
	Instruments in wire baskets inside sterilization container with valve inside lid and bottom.		<b>+</b>	<b>-</b>	<b>+2)</b>	<b>+2)</b>	<b>+2)</b>	<b>+2)</b>
	Instruments in wire baskets inside sterilization container with valve inside lid, no valve in the bottom.		<b>+</b>	<b>-</b>	<b>+2)</b>	<b>+2)</b>	<b>+2)</b>	<b>+2)</b>

Legend: + applicable, – not applicable, 0 applicable under certain conditions

- 1). To avoid service, operation and/or loading errors, see Table 1 for some recommendations regarding the correct or incorrect combination of sterilization containers with standard steam sterilisation procedures.
- 2). Only sterilisation containers from which the condensate evaporation or flow off is guaranteed are suitable.
- 3). The gravity procedure without a drying device does not provide sufficient re-drying.



Standard Model



Safety Model

**GB**

Choose between our standard and our safety model. The safety lid protects against contamination while transport and stocking. Filter- and safety lid can only separate by opening container.

**DE**

Wählen Sie zwischen unserem Standard- und unserem Sicherheitsmodell. Der Sicherheitsdeckel schützt vor Kontamination bei Transport und Lagerung. Filter- und Sicherheitsdeckel lassen sich nur durch Öffnen des Behälters trennen.

**ES**

Elija entre nuestro estándar y nuestro modelo de seguridad. La tapa de seguridad protege contra la contaminación durante el transporte y almacenamiento. El filtro y la tapa de seguridad solo pueden separarse abriendo el recipiente.

**FR**

Choisissez entre notre modèle standard et notre modèle de sécurité. Le couvercle de sécurité protège contre la contamination pendant le transport et le stockage. Le filtre et le couvercle de sécurité ne peuvent se séparer qu'en ouvrant le récipient.

**IT**

Scegli tra il nostro modello standard e il nostro modello di sicurezza. Il coperchio di sicurezza protegge dalla contaminazione durante il trasporto e lo stoccaggio. Il filtro e il coperchio di sicurezza possono essere separati solo aprendo il contenitore.



Lid Perforated, Bottom Non Perforated



Lid Perforated, Bottom Perforated



GB

The locking flaps are movable 180° on a pivot and can be secured in any position.

DE

Die Verschlussklappen sind um 180° auf einem Drehpunkt beweglich und können in jeder Position arretiert werden.

ES

Las aletas de bloqueo se pueden mover 180° sobre un pivote y se pueden asegurar en cualquier posición.

FR

Les volets de verrouillage sont mobiles à 180° sur un pivot et peuvent être bloqués dans n'importe quelle position.

IT

Le alette di chiusura sono mobili di 180° su un perno e possono essere fissate in qualsiasi posizione.



GB

The filter retainer is locked by pressing down in center of the filter. Opening the filter retainer by pressing.

DE

Der Filterhalter wird durch Herunterdrücken in der Mitte des Filters verriegelt. Öffnen des Filterhalters durch Drücken.

ES

El retenedor del filtro se bloquea presionando hacia abajo en el centro del filtro. Apertura del portafiltros presionando.


FR

Le dispositif de retenue du filtre est verrouillé en appuyant au centre du filtre. Ouverture du porte-filtre par pression.

IT

Il fermo del filtro si blocca premendo al centro del filtro. Apertura del portafiltro premendo.


1/1 Container - Standard Model



Non perforated  
Bottom

- 21-444-01
- 21-444-02
- 21-444-03
- 21-444-04
- 21-444-05
- 21-444-06

Weight: 2634 gm



Perforated  
Bottom

- 21-445-01
- 21-445-02
- 21-445-03
- 21-445-04
- 21-445-05
- 21-445-06

Weight: 3156 gm



Outside: 580 x 280 x 100 mm  
 Inside: 550 x 265 x 84 mm

- 21-446-01
- 21-446-02
- 21-446-03
- 21-446-04
- 21-446-05
- 21-446-06

Weight: 3125.5 gm

- 21-447-01
- 21-447-02
- 21-447-03
- 21-447-04
- 21-447-05
- 21-447-06

Weight: 3244 gm



Outside: 580 x 280 x 135 mm  
 Inside: 550 x 265 x 115 mm

- 21-444-01
- 21-444-02
- 21-444-03
- 21-444-04
- 21-444-05
- 21-444-06

Weight: 2634 gm

- 21-445-01
- 21-445-02
- 21-445-03
- 21-445-04
- 21-445-05
- 21-445-06

Weight: 3156 gm



Outside: 580 x 280 x 150 mm  
 Inside: 550 x 265 x 130 mm

- 21-448-01
- 21-448-02
- 21-448-03
- 21-448-04
- 21-448-05
- 21-448-06

Weight: 3219 gm

- 21-449-01
- 21-449-02
- 21-449-03
- 21-449-04
- 21-449-05
- 21-449-06

Weight: 3382 gm



Outside: 580 x 280 x 200 mm  
 Inside: 550 x 265 x 180 mm

- 21-450-01
- 21-450-02
- 21-450-03
- 21-450-04
- 21-450-05
- 21-450-06

Weight: 3429.5 gm

- 21-451-01
- 21-451-02
- 21-451-03
- 21-451-04
- 21-451-05
- 21-451-06

Weight: 3738.5 gm



Outside: 580 x 280 x 260 mm  
 Inside: 550 x 265 x 240 mm


1/1 Container - Safety Model



Non perforated  
Bottom

- 21-454-01
- 21-454-02
- 21-454-03
- 21-454-04
- 21-454-05
- 21-454-06

Weight: 3486 gm



Perforated  
Bottom

- 21-455-01
- 21-455-02
- 21-455-03
- 21-455-04
- 21-455-05
- 21-455-06

Weight: 4007.5 gm



Outside: 580 x 280 x 100 mm  
Inside: 550 x 265 x 84 mm

- 21-456-01
- 21-456-02
- 21-456-03
- 21-456-04
- 21-456-05
- 21-456-06

Weight: 3981.5 gm

- 21-457-01
- 21-457-02
- 21-457-03
- 21-457-04
- 21-457-05
- 21-457-06

Weight: 4095.5 gm



Outside: 580 x 280 x 135 mm  
Inside: 550 x 265 x 115 mm

- 21-458-01
- 21-458-02
- 21-458-03
- 21-458-04
- 21-458-05
- 21-458-06

Weight: 4070.5 gm

- 21-459-01
- 21-459-02
- 21-459-03
- 21-459-04
- 21-459-05
- 21-459-06

Weight: 4233 gm



Outside: 580 x 280 x 150 mm  
Inside: 550 x 265 x 130 mm

- 21-460-01
- 21-460-02
- 21-460-03
- 21-460-04
- 21-460-05
- 21-460-06

Weight: 4280 gm

- 21-461-01
- 21-461-02
- 21-461-03
- 21-461-04
- 21-461-05
- 21-461-06

Weight: 4589 gm



Outside: 580 x 280 x 200 mm  
Inside: 550 x 265 x 180 mm

- 21-462-01
- 21-462-02
- 21-462-03
- 21-462-04
- 21-462-05
- 21-462-06

Weight: 4770.5 gm


- 21-463-01
- 21-463-02
- 21-463-03
- 21-463-04
- 21-463-05
- 21-463-06

Weight: 5140 gm




Outside: 580 x 280 x 260 mm  
Inside: 550 x 265 x 240 mm

1/1 Container - Bottoms and Lids



Non perforated  
Bottom

	Outside	Inside	Weight
21-464-01	580 x 280 x 100 mm	550 x 262 x 84 mm	1923 gm
21-464-02	580 x 280 x 135 mm	550 x 262 x 117 mm	2019 gm
21-464-03	580 x 280 x 150 mm	550 x 262 x 132 mm	2108 gm
21-464-04	580 x 280 x 200 mm	550 x 262 x 182 mm	2317 gm
21-464-05	580 x 280 x 260 mm	550 x 262 x 244 mm	2808 gm



Perforated  
Bottom

	Outside	Inside	Weight
21-465-01	580 x 280 x 100 mm	550 x 262 x 84 mm	2044 gm
21-465-02	580 x 280 x 135 mm	550 x 262 x 117 mm	2132 gm
21-465-03	580 x 280 x 150 mm	550 x 262 x 132 mm	2270 gm
21-465-04	580 x 280 x 200 mm	550 x 262 x 182 mm	2627 gm
21-465-05	580 x 280 x 260 mm	550 x 262 x 244 mm	2851 gm



Perforated lid



- 21-466-01
- 21-466-02
- 21-466-03
- 21-466-04
- 21-466-05
- 21-466-06

Weight: 1111 gm  
 Size: 580 x 280 mm

Safety Lid



- 21-467-01
- 21-467-02
- 21-467-03
- 21-467-04
- 21-467-05
- 21-467-06

Weight: 851.5 gm  
 Size: 580 x 280 mm

1/1 Container - Wire Baskets



Wire Basket

21-548-01	485 x 255 x 30 mm
21-548-02	485 x 255 x 50 mm
21-548-03	485 x 255 x 70 mm
21-548-04	485 x 255 x 100 mm



Wire Basket

21-545-01	540 x 255 x 30 mm
21-545-02	540 x 255 x 50 mm
21-545-03	540 x 255 x 70 mm
21-545-04	540 x 255 x 100 mm



Wire Base

21-611-03	545 x 255 x 30 mm
-----------	-------------------



Drape Retainer

21-612-05	545 x 255 mm
-----------	--------------



Drape

21-610-07	130 x 60 cm
21-610-08	140 x 100 cm
21-610-09	185 x 45 cm
21-610-10	185 x 60 cm


3/4 Container - Safety Model



Non perforated  
Bottom

Grey	21-468-01
Yellow	21-468-02
Red	21-468-03
Blue	21-468-04
Green	21-468-05
Black	21-468-06

Weight: 2268.5 gm



Perforated  
Bottom

Grey	21-469-01
Yellow	21-469-02
Red	21-469-03
Blue	21-469-04
Green	21-469-05
Black	21-469-06

Weight: 2503 gm



Outside: 465 x 280 x 100 mm  
 Inside: 425 x 265 x 84 mm

Grey	21-470-01
Yellow	21-470-02
Red	21-470-03
Blue	21-470-04
Green	21-470-05
Black	21-470-06

Weight: 2549.5 gm

Grey	21-471-01
Yellow	21-471-02
Red	21-471-03
Blue	21-471-04
Green	21-471-05
Black	21-471-06

Weight: 2715 gm



Outside: 465 x 280 x 135 mm  
 Inside: 425 x 265 x 115 mm

Grey	21-472-01
Yellow	21-472-02
Red	21-472-03
Blue	21-472-04
Green	21-472-05
Black	21-472-06

Weight: 2620.5 gm


Grey	21-473-01
Yellow	21-473-02
Red	21-473-03
Blue	21-473-04
Green	21-473-05
Black	21-473-06

Weight: 2718 gm



Outside: 465 x 280 x 150 mm  
 Inside: 425 x 260 x 130 mm


3/4 Container - Safety Model



Non perforated  
Bottom

Grey	21-474-01
Yellow	21-474-02
Red	21-474-03
Blue	21-474-04
Green	21-474-05
Black	21-474-06

Weight: 2980 gm



Perforated  
Bottom

Grey	21-475-01
Yellow	21-475-02
Red	21-475-03
Blue	21-475-04
Green	21-475-05
Black	21-475-06

Weight: 3213.5 gm



Outside: 465 x 280 x 100 mm  
Inside: 425 x 265 x 84 mm

Grey	21-476-01
Yellow	21-476-02
Red	21-476-03
Blue	21-476-04
Green	21-476-05
Black	21-476-06

Weight: 3260 gm

Grey	21-477-01
Yellow	21-477-02
Red	21-477-03
Blue	21-477-04
Green	21-477-05
Black	21-477-06

Weight: 2425 gm



Outside: 465 x 280 x 135 mm  
Inside: 425 x 265 x 115 mm

Grey	21-478-01
Yellow	21-478-02
Red	21-478-03
Blue	21-478-04
Green	21-478-05
Black	21-478-06

Weight: 2975 gm


Grey	21-479-01
Yellow	21-479-02
Red	21-479-03
Blue	21-479-04
Green	21-479-05
Black	21-479-06

Weight: 3561 gm




Outside: 465 x 280 x 150 mm  
Inside: 425 x 260 x 130 mm

3/4 Container - Bottoms and Lids



Non perforated  
 Bottom

	Outside	Inside	Weight
21-480-01	465 x 280 x 100 mm	465 x 262 x 84 mm	1396 gm
21-480-02	465 x 280 x 135 mm	465 x 262 x 117 mm	1676 gm
21-480-03	465 x 280 x 150 mm	465 x 262 x 132 mm	1792 gm



Perforated  
 Bottom

	Outside	Inside	Weight
21-481-01	465 x 280 x 100 mm	425 x 262 x 84 mm	2044 gm
21-481-02	465 x 280 x 135 mm	425 x 262 x 117 mm	2132 gm
21-481-03	465 x 280 x 150 mm	425 x 262 x 132 mm	2270 gm



Perforated lid



- 21-482-01
- 21-482-02
- 21-482-03
- 21-482-04
- 21-482-05
- 21-482-06

Weight: 873 gm  
 Size: 465 x 280 mm

Safety Lid



- 21-483-01
- 21-483-02
- 21-483-03
- 21-483-04
- 21-483-05
- 21-483-06

Weight: 710 gm  
 Size: 465 x 280 mm

### 3/4 Container - Wire Baskets



Wire Basket

21-551-01	405 x 255 x 30 mm
21-551-02	405 x 255 x 50 mm
21-551-03	405 x 255 x 70 mm
21-551-04	405 x 255 x 100 mm



Wire Base

21-611-02	405 x 255 x 30 mm
-----------	-------------------



Drape Retainer


21-612-04	405 x 255 mm
-----------	--------------




Drape

21-610-08	140 x 100 cm
-----------	--------------

1/2 Container - Standard Model



Non perforated  
Bottom



Perforated  
Bottom

- 21-484-01
- 21-484-02
- 21-484-03
- 21-484-04
- 21-484-05
- 21-484-06

Weight: 1744.5 gm

- 21-485-01
- 21-485-02
- 21-485-03
- 21-485-04
- 21-485-05
- 21-485-06

Weight: 1890.5 gm



Outside: 285 x 280 x 100 mm  
 Inside: 262 x 260 x 84 mm

- 21-486-01
- 21-486-02
- 21-486-03
- 21-486-04
- 21-486-05
- 21-486-06

Weight: 1947 gm

- 21-487-01
- 21-487-02
- 21-487-03
- 21-487-04
- 21-487-05
- 21-487-06

Weight: 2125 gm



Outside: 285 x 280 x 135 mm  
 Inside: 262 x 260 x 115 mm

- 21-488-01
- 21-488-02
- 21-488-03
- 21-488-04
- 21-488-05
- 21-488-06

Weight: 1981 gm

- 21-489-01
- 21-489-02
- 21-489-03
- 21-489-04
- 21-489-05
- 21-489-06

Weight: 2462.5 gm



Outside: 285 x 280 x 150 mm  
 Inside: 262 x 260 x 132 mm

- 21-490-01
- 21-490-02
- 21-490-03
- 21-490-04
- 21-490-05
- 21-490-06

Weight: 2101 gm

- 21-491-01
- 21-491-02
- 21-491-03
- 21-491-04
- 21-491-05
- 21-491-06

Weight: 2513.5 gm



Outside: 285 x 280 x 200 mm  
 Inside: 262 x 260 x 180 mm

- 21-492-01
- 21-492-02
- 21-492-03
- 21-492-04
- 21-492-05
- 21-492-06

Weight: 2483.5 gm

- 21-493-01
- 21-493-02
- 21-493-03
- 21-493-04
- 21-493-05
- 21-493-06

Weight: 2657.5 gm




Outside: 285 x 280 x 260 mm  
 Inside: 262 x 260 x 238 mm

1/2 Container - Safety Model



Non perforated  
Bottom



Perforated  
Bottom

Grey	21-494-01
Yellow	21-494-02
Red	21-494-03
Blue	21-494-04
Green	21-494-05
Black	21-494-06

Weight: 2162.5 gm

Grey	21-495-01
Yellow	21-495-02
Red	21-495-03
Blue	21-495-04
Green	21-495-05
Black	21-495-06

Weight: 2309.5 gm



Outside: 285 x 280 x 100 mm  
Inside: 262 x 260 x 84 mm

Grey	21-496-01
Yellow	21-496-02
Red	21-496-03
Blue	21-496-04
Green	21-496-05
Black	21-496-06

Weight: 2351 gm

Grey	21-497-01
Yellow	21-497-02
Red	21-497-03
Blue	21-497-04
Green	21-497-05
Black	21-497-06

Weight: 2466 gm



Outside: 285 x 280 x 135 mm  
Inside: 262 x 260 x 116 mm

Grey	21-498-01
Yellow	21-498-02
Red	21-498-03
Blue	21-498-04
Green	21-498-05
Black	21-498-06

Weight: 2394 gm

Grey	21-499-01
Yellow	21-499-02
Red	21-499-03
Blue	21-499-04
Green	21-499-05
Black	21-499-06

Weight: 2583 gm



Outside: 285 x 280 x 150 mm  
Inside: 262 x 260 x 132 mm

Grey	21-500-01
Yellow	21-500-02
Red	21-500-03
Blue	21-500-04
Green	21-500-05
Black	21-500-06

Weight: 2560 gm

Grey	21-501-01
Yellow	21-501-02
Red	21-501-03
Blue	21-501-04
Green	21-501-05
Black	21-501-06

Weight: 2785.5 gm



Outside: 285 x 280 x 200 mm  
Inside: 262 x 260 x 182 mm

Grey	21-502-01
Yellow	21-502-02
Red	21-502-03
Blue	21-502-04
Green	21-502-05
Black	21-502-06

Weight: 2930 gm

Grey	21-503-01
Yellow	21-503-02
Red	21-503-03
Blue	21-503-04
Green	21-503-05
Black	21-503-06

Weight: 3069.5 gm



Outside: 285 x 280 x 260 mm  
Inside: 262 x 260 x 238 mm

1/2 Container - Bottoms and Lids



Non perforated  
Bottom

	Outside	Inside	Weight
21-504-01	285 x 280 x 100 mm	262 x 262 x 82 mm	1096 gm
21-504-02	285 x 280 x 135 mm	262 x 262 x 120 mm	1284 gm
21-504-03	285 x 280 x 150 mm	262 x 262 x 135 mm	1327 gm
21-504-04	285 x 280 x 200 mm	262 x 262 x 185 mm	1447 gm
21-504-05	285 x 280 x 260 mm	262 x 262 x 240 mm	1817 gm



Perforated  
Bottom

	Outside	Inside	Weight
21-505-01	285 x 280 x 100 mm	262 x 262 x 82 mm	1242 gm
21-505-02	285 x 280 x 135 mm	262 x 262 x 120 mm	1398 gm
21-505-03	285 x 280 x 150 mm	262 x 262 x 135 mm	1471 gm
21-505-04	285 x 280 x 200 mm	262 x 262 x 185 mm	1672 gm
21-505-05	285 x 280 x 260 mm	262 x 262 x 240 mm	2000 gm



Perforated lid



- 21-506-01
- 21-506-02
- 21-506-03
- 21-506-04
- 21-506-05
- 21-506-06

Weight: 254 gm  
 Size: 285 x 280 mm

Safety Lid



- 21-507-01
- 21-507-02
- 21-507-03
- 21-507-04
- 21-507-05
- 21-507-06

Weight: 458 gm  
 Size: 285 x 280 mm

1/2 Container - Wire Baskets



Wire Basket

21-554-01	255 x 245 x 30 mm
21-554-02	255 x 245 x 50 mm
21-554-03	255 x 245 x 70 mm
21-554-04	255 x 245 x 100 mm



Wire Base

21-611-01	255 x 245 x 30 mm
-----------	-------------------



Drape Retainer


21-612-03	255 x 245 mm
-----------	--------------




Drape

21-610-04	90 x 60 cm
21-610-05	90 x 80 cm
21-610-06	130 x 45 cm
21-610-07	130 x 60 cm

Flat Container - Containers and Bottoms



Non perforated  
Bottom



Perforated  
Bottom

- 21-508-01
- 21-508-02
- 21-508-03
- 21-508-04
- 21-508-05
- 21-508-06

Weight: 1285 gm



Outside: 285 x 280 x 55 mm  
 Inside: 270 x 270 x 33 mm

- 21-509-01
- 21-509-02
- 21-509-03
- 21-509-04
- 21-509-05
- 21-509-06

Weight: 1377 gm

- 21-510-01
- 21-510-02
- 21-510-03
- 21-510-04
- 21-510-05
- 21-510-06

Weight: 1529 gm



Outside: 600 x 400 x 180 mm  
 Inside: 550 x 350 x 160 mm



Non perforated  
Bottom

	Outside	Inside	Weight
21-511-01	285 x 280 x 55 mm	270 x 270 x 33 mm	611 gm
21-511-02	285 x 280 x 85 mm	270 x 270 x 64 mm	703 gm




Perforated  
Bottom

	Outside	Inside	Weight
21-512-01	285 x 280 x 85 mm	270 x 270 x 64 mm	857 gm



## Flat Container - Lids and accessories



Wire Basket

21-554-01	255 x 245 x 30 mm
21-554-02	255 x 245 x 50 mm




Drape Retainer

21-612-03	255 x 245 mm
-----------	--------------



Wire Base

21-611-01	255 x 245 x 30 mm
-----------	-------------------



Silicone Mat


21-628-02	245 x 245 mm
-----------	--------------



Drape

21-610-04	90 x 60 cm
21-610-05	90 x 80 cm


## Mini Container



Non perforated  
Bottom

Grey	21-514-01
Yellow	21-514-02
Red	21-514-03
Blue	21-514-04
Green	21-514-05
Black	21-514-06

Weight: 510.5 gm



Perforated  
Bottom

Grey	21-515-01
Yellow	21-515-02
Red	21-515-03
Blue	21-515-04
Green	21-515-05
Black	21-515-06

Weight: 589 gm



Outside: 300 x 140 x 40 mm  
 Inside: 280 x 130 x 27 mm

Grey	21-516-01
Yellow	21-516-02
Red	21-516-03
Blue	21-516-04
Green	21-516-05
Black	21-516-06

Weight: 567 gm

Grey	21-517-01
Yellow	21-517-02
Red	21-517-03
Blue	21-517-04
Green	21-517-05
Black	21-517-06

Weight: 630 gm



Outside: 300 x 140 x 70 mm  
 Inside: 280 x 130 x 55 mm

Grey	21-518-01
Yellow	21-518-02
Red	21-518-03
Blue	21-518-04
Green	21-518-05
Black	21-518-06

Weight: 625.5 gm

Grey	21-519-01
Yellow	21-519-02
Red	21-519-03
Blue	21-519-04
Green	21-519-05
Black	21-519-06

Weight: 625 gm



Outside: 300 x 140 x 100 mm  
 Inside: 280 x 130 x 84 mm



Non perforated  
Bottom

	Outside	Inside	Weight
21-520-01	300 x 140 x 40 mm	280 x 130 x 27 mm	190 gm
21-520-02	300 x 140 x 70 mm	280 x 130 x 55 mm	241 gm
21-520-03	300 x 140 x 100 mm	280 x 130 x 84 mm	305 gm




Perforated  
Bottom

	Outside	Inside	Weight
21-521-01	300 x 140 x 40 mm	280 x 130 x 27 mm	258 gm
21-521-02	300 x 140 x 70 mm	280 x 130 x 55 mm	304 gm
21-521-03	300 x 140 x 100 mm	280 x 130 x 84 mm	326 gm

## Mini Container - Bottoms and Lids



Perforated lid



Grey	21-522-01
Yellow	21-522-02
Orange	21-522-03
Red	21-522-04
Blue	21-522-05
Green	21-522-06
Black	21-522-06

Weight: 325.5 gm

## Mini Container - Accessories



Drape Retainer

21-612-01 280 x 125 mm



Silicone Mat


21-628-01 270 x 125 mm



Drape

21-610-01 60 x 40 cm

## Dental Containers




Perforated  
Bottom

21-523-01

Weight: 783 gm



Outside: 310 x 190 x 40 mm  
Inside: 280 x 182 x 35 mm



Perforated  
Bottom

21-524-01

Weight: 862 gm



Outside: 310 x 190 x 65 mm  
Inside: 280 x 182 x 55 mm



Perforated  
Bottom


21-525-01

Weight: 1022 gm



Outside: 310 x 190 x 130 mm  
Inside: 280 x 182 x 124 mm

Dental Container - Bottom, Lid and Accessories



Perforated Bottom

	Outside	Inside	Weight
21-526-01	310 x 190 x 40 mm	280 x 182 x 35 mm	341 gm
21-526-02	310 x 190 x 65 mm	280 x 182 x 55 mm	415 gm
21-526-03	310 x 190 x 130 mm	280 x 182 x 124 mm	570 gm



Perforated lid



21-527-01

Weight: 417 gm  
 Size: 310 x 190 mm



Drape Retainer

21-612-02 310 x 190 mm



Drape

21-610-01	60 x 40 cm
21-610-02	65 x 65 cm
21-610-03	85 x 85 cm

## Implant Containers



Lid perforated  
Bottom perforated (no filter)

Grey	21-528-01
Yellow	21-528-02
Red	21-528-03
Blue	21-528-04
Green	21-528-05
Black	21-528-06

Weight: 654 gm



Outside: 500 x 155 x 75 mm  
Inside: 495 x 155 x 75 mm



Lid perforated  
Bottom perforated (no filter)

Grey	21-529-01
Yellow	21-529-02
Red	21-529-03
Blue	21-529-04
Green	21-529-05
Black	21-529-06

Weight: 368 gm



Outside: 500 x 155 x 75 mm  
Inside: 495 x 155 x 75 mm

## Optik Containers



Lid perforated  
Bottom perforated (no filter)

21-530-01

Weight: 302 gm



Outside: 450 x 70 x 70 mm  
Inside: 450 x 70 x 30 mm



Lid perforated  
Bottom perforated (no filter)

21-531-01

Weight: 152 gm



Outside: 250 x 60 x 50 mm  
Inside: 250 x 55 x 50 mm

## Removal Containers



Non perforated  
Lid and Bottom

Grey	21-532-01
Yellow	21-532-02
Red	21-532-03
Blue	21-532-04
Green	21-532-05
Black	21-532-06

Weight: 1677 gm



Outside: 285 x 280 x 100 mm  
Inside: 262 x 260 x 82 mm




Non perforated  
Lid and Bottom

Grey	21-533-01
Yellow	21-533-02
Red	21-533-03
Blue	21-533-04
Green	21-533-05
Black	21-533-06

Weight: 1929 gm



Outside: 285 x 280 x 150 mm  
Inside: 262 x 260 x 132 mm



Non perforated  
Lid and Bottom


Grey	21-534-01
Yellow	21-534-02
Red	21-534-03
Blue	21-534-04
Green	21-534-05
Black	21-534-06

Weight: 2883 gm




Outside: 580 x 280 x 150 mm  
Inside: 548 x 258 x 130 mm

Removal Container - Bottoms and Lids



Non perforated  
Bottom

	Outside	Weight
21-535-01	580 x 280 x 150 mm	2108 gm



Non perforated  
Bottom

	Outside	Weight
21-536-01	285 x 280 x 150 mm	1096 gm
21-636-02	285 x 280 x 150 mm	1327 gm



Non-Perforated lid



21-537-01
21-537-02
21-537-03
21-537-04
21-537-05
21-537-06

Weight: 458 gm  
 Size: 285 x 280 mm

Non-Perforated lid



21-538-01
21-538-02
21-538-03
21-538-04
21-538-05
21-538-06

Weight: 875 gm  
 Size: 580 x 280 mm

## Brush Container - Stainless Steel and Aluminium



21-539-10 For 10 brushes

21-539-20 For 20 brushes

Stainless steel

■	21-540-01
■	21-540-02
■	21-540-03
■	21-540-04
■	21-540-05

Aluminium  
For 20 brushes  
Size: 240 x 275 x 70 mm



#### Single Tray Containers

21-541-01	200 x 140 x 50 mm
21-541-02	280 x 195 x 50 mm
21-541-03	265 x 255 x 50 mm
21-541-04	380 x 260 x 50 mm

For micro instruments



#### Double Tray Containers

21-542-01	280 x 195 x 50 mm
21-542-02	265 x 255 x 50 mm
21-542-03	380 x 260 x 50 mm

With 2 take-out perforated trays

## Baskets



21-543-01	540 x 255 x 30 mm
21-543-02	540 x 255 x 50 mm
21-543-03	540 x 255 x 70 mm
21-543-04	540 x 255 x 100 mm

Wire basket with lid

21-544-00

Lid only



21-545-01	540 x 255 x 30 mm
21-545-02	540 x 255 x 50 mm
21-545-03	540 x 255 x 70 mm
21-545-04	540 x 255 x 100 mm

Wire basket



21-546-01	485 x 255 x 30 mm
21-546-02	485 x 255 x 50 mm
21-546-03	485 x 255 x 70 mm
21-546-04	485 x 255 x 100 mm

Wire basket with lid

21-547-00

Lid only



21-548-01	485 x 255 x 30 mm
21-548-02	485 x 255 x 50 mm
21-548-03	485 x 255 x 70 mm
21-548-04	485 x 255 x 100 mm

Wire basket

## Baskets



21-549-01	405 x 255 x 30 mm
21-549-02	405 x 255 x 50 mm
21-549-03	405 x 255 x 70 mm
21-549-04	405 x 255 x 100 mm

Wire basket with lid

21-550-00

Lid only



21-551-01	405 x 255 x 30 mm
21-551-02	405 x 255 x 50 mm
21-551-03	405 x 255 x 70 mm
21-551-04	405 x 255 x 100 mm

Wire basket



21-552-01	255 x 245 x 30 mm
21-552-02	255 x 245 x 50 mm
21-552-03	255 x 245 x 70 mm
21-552-04	255 x 245 x 100 mm

Wire basket with lid

21-553-00

Lid only



21-554-01	255 x 245 x 30 mm
21-554-02	255 x 245 x 50 mm
21-554-03	255 x 245 x 70 mm
21-554-04	255 x 245 x 100 mm

Wire basket

## Baskets



21-555-01 540 x 255 x 50 mm

21-555-02 540 x 255 x 70 mm

21-555-03 540 x 255 x 100 mm

Wire basket with lid

21-556-00

Lid only

21-557-01 485 x 255 x 50 mm

21-557-02 485 x 255 x 70 mm

21-557-03 485 x 255 x 100 mm

Wire basket with lid

21-558-00

Lid only



21-561-01 405 x 255 x 50 mm

21-561-02 405 x 255 x 70 mm

21-561-03 405 x 255 x 100 mm

Wire basket with lid

21-562-00

Lid only



21-564-01 255 x 245 x 50 mm


21-564-02 255 x 245 x 70 mm

21-564-03 255 x 245 x 100 mm

Wire basket with lid

21-565-00

Lid only



21-559-01 540 x 255 x 50 mm

21-559-02 540 x 255 x 70 mm

21-559-03 540 x 255 x 100 mm

Wire basket

21-560-01 485 x 255 x 50 mm

21-560-02 485 x 255 x 70 mm

21-560-03 485 x 255 x 100 mm

Wire basket



21-563-01 405 x 255 x 50 mm

21-563-02 405 x 255 x 70 mm

21-563-03 405 x 255 x 100 mm

Wire basket




21-566-01 255 x 245 x 50 mm

21-566-02 255 x 245 x 70 mm

21-566-03 255 x 245 x 100 mm


Wire basket

## Baskets



21-567-01 290 x 80 x 52 mm  
 21-567-02 460 x 80 x 52 mm  
 21-567-03 670 x 80 x 52 mm

Endoscope basket



21-568-00 640 x 150 x 77 mm

Endoscope basket



21-569-00 274 x 172 x 43 mm

Wire basket



21-570-00 267 x 125 x 50 mm

Wire basket



21-571-00 274 x 172 x 43 mm


Wire basket



21-572-00 267 x 125 x 50 mm

Wire basket

## Wire Baskets



21-573-01	30 x 30 x 30 mm
21-573-02	40 x 40 x 20 mm
21-573-03	40 x 40 x 40 mm
21-573-04	60 x 60 x 20 mm
21-573-05	60 x 60 x 40 mm
21-573-06	60 x 60 x 60 mm
21-573-07	80 x 80 x 20 mm
21-573-08	80 x 40 x 40 mm
21-573-09	80 x 80 x 40 mm
21-573-10	80 x 80 x 60 mm
21-573-11	80 x 80 x 80 mm
21-573-12	100 x 100 x 10 mm
21-573-13	120 x 120 x 30 mm
21-573-14	120 x 120 x 60 mm
21-573-15	120 x 120 x 80 mm
21-573-16	240 x 240 x 30 mm
21-573-17	240 x 240 x 60 mm
21-573-18	240 x 240 x 80 mm
21-573-19	240 x 240 x 100 mm

Small fine mesh basket with lid  
Mesh distance, 2 x 2 mm



21-574-01 600 x 400 x 150 mm

21-574-02 600 x 400 x 300 mm

Laundry basket



21-575-01 600 x 400 x 150 mm

21-575-02 600 x 400 x 300 mm

Laundry basket

## Baskets



Size: 1/1

Box and lid	Box only	Lid only	
21-576-01	21-577-01	21-578-01	540 x 253 x 30 mm
21-576-02	21-577-02	21-578-02	540 x 253 x 50 mm
21-576-03	21-577-03	21-578-03	540 x 253 x 70 mm
21-576-04	21-577-04	21-578-04	540 x 253 x 100 mm



Size: 1/1

Box and lid	Box only	Lid only	
21-579-01	21-580-01	21-581-01	485 x 253 x 30 mm
21-579-02	21-580-02	21-581-02	485 x 253 x 50 mm
21-579-03	21-580-03	21-581-03	485 x 253 x 70 mm
21-579-04	21-580-04	21-581-04	485 x 253 x 100 mm



Size: 3/4

Box and lid	Box only	Lid only	
21-582-01	21-583-01	21-584-01	405 x 253 x 30 mm
21-582-02	21-583-02	21-584-02	405 x 253 x 50 mm
21-582-03	21-583-03	21-584-03	405 x 253 x 70 mm
21-582-04	21-583-04	21-584-04	405 x 253 x 100 mm



Size: 1/2

Box and lid	Box only	Lid only	
21-585-01	21-586-01	21-587-01	244 x 253 x 30 mm
21-585-02	21-586-02	21-587-02	244 x 253 x 50 mm
21-585-03	21-586-03	21-587-03	244 x 253 x 70 mm
21-585-04	21-586-04	21-587-04	244 x 253 x 100 mm



Size: Dental

Box and lid	Box only	Lid only	
21-588-01	21-589-01	21-590-01	275 x 175 x 40 mm
21-588-02	21-589-02	21-590-02	275 x 175 x 95 mm



Size: Small

Box and lid	Box only	Lid only	
21-591-01	21-592-01	21-593-01	265 x 125 x 50 mm



21-594-14  
 Stand  
 14 mm height




21-595-20  
 Stand  
 20 mm height

## Filters



**21-596-00** 210x95 mm

Long term textile filter for mini and dental container for 50 sterilization cycles.



**21-597-00** 235 x 118 mm

Long term PTFE filter for mini and dental container for 1200 sterilization cycles.



**21-598-00** 235 x 118 mm

Long term PTFE filter for 1/1, 1/2, 3/4, mini and dental container.



**21-599-00** Ø 190 mm

Long term textil filter for 1/1, 1/2, 3/4 container for 50 sterilization cycles.



**21-600-00** 235 x 118 mm

Long term PTFE filter for 1/1, 1/2, 3/4 container for 1200 sterilization cycles.



**21-601-00** 235 x 118 mm

Disposable filter for 1/1, 1/2, 3/4 container.

Indicators, Labels and Drapes



21-602-00

Paper label with indicator  
 For 1/1, 1/2, 3/4 containers  
 Pack of 100



21-602-00

Paper label with indicator  
 For mini and dental containers  
 Pack of 100



21-602-00

Heavy duty clip  
 Serrated jaws



21-605-01 For paper filter  
 21-605-02 For textile filter

Filter tray



21-606-00

Filter dispenser




21-607-00

Plastic security seal  
 Pack of 100



21-607-00

Dispenser for  
 identification labels



Coloured identification labels.  
 Printed to your specification,  
 16 letters max.



Drape

21-610-01	60 x 40 cm
21-610-02	65 x 65 cm
21-610-03	85 x 85 cm
21-610-04	90 x 60 cm
21-610-05	90 x 80 cm
21-610-06	130 x 45 cm
21-610-07	130 x 60 cm
21-610-08	140 x 100 cm
21-610-09	185 x 45 cm
21-610-10	185 x 60 cm

## Wire Base, Drape Retainer, Metal Labels and Fixation Clamps



Wire Base

21-611-01	255 x 245 x 30 mm
21-611-02	405 x 255 x 30 mm
21-611-03	545 x 255 x 30 mm



Drape Retainer

21-612-01	280 x 125 mm
21-612-02	310 x 190 mm
21-612-03	255 x 245 mm
21-612-04	405 x 255 mm
21-612-05	545 x 255 mm



Metal Label for Wire Basket

21-613-00	545 x 255 x 30 mm
-----------	-------------------

Unlabeled



Metal Label for Wire Basket

21-614-00	545 x 255 mm
-----------	--------------

Labeled



21-615-01

24 x 28 x 25 mm  
Ø 7 - 12 mm



21-615-02

24 x 28 x 29 mm  
Ø 12 - 16 mm



21-615-03

24 x 28 x 31 mm  
Ø 16 - 26 mm



21-615-04

24 x 28 x 40 mm  
Ø 26 - 36 mm



21-615-05

24 x 28 x 44 mm  
Ø 36 - 45 mm




21-616-01

50 x 20 x 12 mm



21-616-02

65 x 20 x 20 mm



## Tray Tags



	21-617-01
	21-617-02
	21-617-03
	21-617-04
	21-617-05

	21-617-01
	21-617-02
	21-617-03
	21-617-04
	21-617-05

Nylon  
 Size: 105x40mm  
 Pack of 100

21-618-00

Stainless steel  
 Attachment ring  
 Ø 30 mm  
 Pack of 100

## Components




21-619-00

Holdong pin for  
partition sheets




21-620-00

Holdong pin for  
wire baskets




21-621-01	50 x 20 mm
21-621-02	130 x 20 mm
21-621-03	225 x 20 mm
21-621-04	460 x 20 mm

Sheets



21-622-01	230 mm
21-622-02	470 mm

Sliding rails with securing screws



21-623-01	230 mm
21-623-02	470 mm

Limiting bar with securing screws



21-624-00

Divider, to be used  
with 21-622-01 and 21-622-02



21-625-00


Spiral holding



21-626-00


Double spiral holding

## Silicone Mats




Silicone Mat Fig. 1

21-627-01	270 x 125 mm
21-627-02	245 x 245 mm
21-627-03	220 x 230 mm
21-627-04	380 x 230 mm
21-627-05	440 x 230 mm
21-627-06	520 x 230 mm




Silicone Mat Fig. 2

21-628-01	270 x 125 mm
21-628-02	245 x 245 mm
21-628-03	220 x 230 mm
21-628-04	380 x 230 mm
21-628-05	440 x 230 mm
21-628-06	520 x 230 mm



Silicone Mat Fig. 3

21-629-00	545 x 255 mm
-----------	--------------



Silicone Mat Fig. 4


21-630-01	220 x 125 mm
21-630-02	380 x 230 mm
21-630-03	440 x 230 mm
21-630-04	520 x 230 mm

## Silicone Fixation System




Length Height

21-631-01	97 mm	41 mm
21-631-02	125 mm	41 mm
21-631-03	152 mm	41 mm
21-631-04	250 mm	41 mm





Length Height

21-632-01	90 mm	41 mm
21-632-02	120 mm	41 mm
21-632-03	150 mm	41 mm
21-632-04	240 mm	41 mm




Length Height

21-633-01	87 mm	41 mm
21-633-02	119 mm	41 mm
21-633-03	150 mm	41 mm
21-633-04	238 mm	41 mm



Length Height


21-634-01	77 mm	41 mm
21-634-02	112 mm	41 mm
21-634-03	158 mm	41 mm
21-634-04	238 mm	41 mm



Length Height

21-635-00	178 mm	20 mm
-----------	--------	-------


Instruments holder from Ø 3 to 8 mm



Length Height

21-636-00	73 mm	20 mm
-----------	-------	-------


Instruments holder from Ø 5



Length Height

21-637-00	73 mm	20 mm
-----------	-------	-------


Instruments holder from Ø 5



Length Height

21-638-00	178 mm	30 mm
-----------	--------	-------

Instruments holder from Ø 3 to 8 mm




## Silicone Fixation System





Length Height

21-639-00 121 mm 18 mm



Length Height


21-640-00 123 mm 18 mm



Length Height

21-641-00 178 mm 20 mm



Instruments holder from Ø 3 mm



Length Height



21-642-00 123 mm 18 mm

Instruments holder from Ø 5 mm



Length Height



21-643-00 123 mm 18 mm



Length Height

21-644-00 23.5 mm 21 mm



Ø 8 mm



Length Height

21-645-00 23.5 mm 21 mm



Ø 5 mm



Length Height

21-646-00 23.5 mm 21 mm


Ø 13 mm




Length Height

21-647-00 23.5 mm 21 mm

Ø 10.5 mm




## Silicone Fixation System




	Length	Height
21-648-00	73 mm	30 mm

Endoscope holder, double  
one side  $\varnothing$  15 mm open  
one side from  $\varnothing$  3 to  $\varnothing$  15 mm




	Length	Height
21-649-00	72 mm	30 mm

Endoscope holder, double  
one side  $\varnothing$  15 mm open  
one side from  $\varnothing$  3 to  $\varnothing$  15 mm




	Length	Height
21-650-00	72 mm	30 mm

Instruments holder  
from  $\varnothing$  3 to  $\varnothing$  15 mm




	Length	Height
21-651-00	72 mm	30 mm

Instruments holder  
from  $\varnothing$  3 to  $\varnothing$  8 mm



	Length	Height
21-652-00	38 mm	20 mm


Endoscope holder, double  
one side 12 mm open  
one side from  $\varnothing$  3 to  $\varnothing$  8 mm



	Length	Height
21-638-00	72 mm	30 mm


Endoscope holder, double  
one side 12 mm open  
one side from  $\varnothing$  3 to  $\varnothing$  8 mm

## Silicone Fixation System




	Length	Height
21-654-00	54 mm	20 mm

Endoscope holder  
 Ø 4 mm closed




	Length	Height
21-655-00	38 mm	29 mm

Endoscope holder  
 Ø 6 mm closed




	Length	Height
21-656-00	38 mm	29 mm

Endoscope holder  
 Ø 8 mm closed




	Length	Height
21-657-00	38 mm	29 mm

Endoscope holder  
 Ø 12 mm closed



	Length	Height
21-658-00	38 mm	29 mm

Endoscope holder  
 Ø 12 mm open



	Length	Height
21-659-00	38 mm	29 mm

Endoscope holder  
 Ø 15 mm open